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SensiTemp® Resin Temporary Cement

Indications For Use

SensiTemp® Resin Temporary Cement is indicated for the provisional (reversible) cementation of temporary dental crowns and bridges.

Contraindications

SensiTemp® Resin Temporary Cement is contraindicated for the permanent cementation of any fixed dental restoration. Due to its high compressive strength, *SensiTemp®* Resin is contraindicated for provisional cementation of temporary intracoronal (inlays and onlays) restorations.

Warnings

SensiTemp® Resin should not be used with patients who have known sensitivities to acrylics, bisacrylics, composite resins, or methacrylates.

Precautions

Care should be taken when using *SensiTemp®* Resin Temporary Cement over or next to new composite resins or resin core build-up materials in order to prevent inadvertent bonding of the cement to the build-up by using a non-reactive lubricant over any composite or resin build-up materials. See Directions for Use for complete instructions. Do not use if product is past expiration date.

Adverse Reactions

SensiTemp® Resin Temporary Cement contains methacrylates. Do not use this product with patients who have known or suspected sensitivities or allergies to methacrylates, acrylics, bisacrylics or composite resins.

Directions for use

5 ml Hand Syringe - First Time Use - Operating the Syringe

- Remove the protective brown cap and discard.
- Gently push on the plunger and express a small amount of cement, ensuring material is flowing equally from both syringe orifices. Carefully wipe clean.
- Assemble a mixing tip to the syringe and turn clock-wise 1/4 turn to lock in place.
- Push on the plunger to express cement through the mixing tip. Discard the first pea-sized amount.

The syringe is now ready to dispense cement directly into a crown. After use, leave the mixing tip in place until the next use. Subsequent use should not require bleeding of the syringe, however for best results we recommend wiping any material from the syringe orifice prior to placing a new mixing tip.

25ml Cartridge – First Time Use – Operating the Cartridge

- Load the cartridge into the dispensing gun ensuring the cartridge is fully seated in the dispenser and the locking latch is fully engaged.
- Remove the blue-colored protective cap from the cartridge and discard.
- Push the black plunger on the dispensing gun so that is seats snugly against the cartridge pistons.
- Squeeze the trigger on the dispenser to bleed a small amount of material from the cartridge, ensuring that both sides are flowing evenly. Using a tissue, wipe the cartridge orifices and discard the bled material.
- Affix the appropriate mixing tip (light blue 3.2 mm) onto the cartridge by aligning the notch on the tip with the notch on the cartridge. Ensure the tip is fully seated, then rotate the blue hub on the tip ¼ turn to lock it into place.
- Squeeze the trigger on the dispenser to extrude the temporary cement. Discard the first pea-sized amount of material from the mixing tip.

The cartridge is now ready to dispense cement directly into a crown. After use, leave the mixing tip in place until the next use. Subsequent use should not require bleeding of the cartridge; however, for best results we recommend wiping any material from the cartridge orifice prior to placing a new mixing tip.

Cementing A Restoration

- Ensure that the temporary restoration fits well, is polished (or coated with GlossCote® Varnish), and the interior is clean and dry.
- Isolate and dry the prepared tooth. *SensiTemp®* Resin may bond to resin core materials. If the tooth has a composite resin build-up or core, very lightly lubricate the preparation with a lubricant. We recommend a thin layer of MQ Lubricant or petrolatum to prevent adhesion to the build-up material.
- Assemble a new mixing tip onto the *SensiTemp®* Resin syringe, and express cement into the crown. Keep the tip buried in the cement flow to prevent air bubbles. Ensure the interior of the crown is coated with cement while not overfilling. This can be achieved by starting with a shallow fill in the provisional crown and using the tip of the syringe/cartridge to pull *SensiTemp®* Resin up to cover the margins. Begin keeping time.
- Seat the crown and have the patient firmly hold in place by closing on a cotton roll.
- At 90-120 seconds, tease the excess with an instrument (probe, explorer or scaler). The chemistry of *SensiTemp®* Resin is such that it will be in a “rubber” phase at

this point, allowing easy removal of the excess cement. Make sure to pass the probe interproximally to remove any interproximal cement, following with dental floss. Do not allow the cement to harden past this “rubber” phase before starting removal of excess cement or removal will be more difficult due to the compressive strength of the product.

- Final cement set is in 4 to 5 minutes.
- If the crown needs to be removed and re-cemented, the film thickness of *SensiTemp®* Resin is thin enough that the initial layer may not need to be removed. If for any reason the cement needs to be removed before re-cementation, use a small, round bur to remove the cement from the crown interior.

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